



Stability constants and thermodynamic properties of complexation of Ibuprofen with Zr^{4+} PH metrically

G. D. Rawate*

*Department of Chemistry, Shri R.R.Lahoti Science College, Morshi, Dist-Amravati, Maharashtra, India

Abstract : pHmetric studies of binary complex of Zn^{4+} with ibuprofen has been carried out in ethanol at three temperatures (25 ± 0.1 , 30 ± 0.1 and 35 ± 0.1 °C) and at an ionic strength of 0.1 mol L^{-1} (KNO_3). The thermodynamic parameters ΔG , ΔH and ΔS are calculated by known equation. pHmetric measurement of hydrogen ion concentration may be employed when the degree of complex formation is sensitive to the hydrogen ion concentration thus the degree of complex formation undergoes increase/ decrease with change in pH. The method of Calvin and Bjerrum as adopted by Irving and Rossotti has been employed to determine log K1 values. The thermodynamic parameters ΔG , ΔH and ΔS are calculated.

Keywords : pH metry, Stability constant, Thermodynamic parameters, Free energy, Enthalpy, Entropy.

Introduction:

Ibuprofen is used primarily for fever. Ibuprofen, from isobutylphenylpropanoic acid, is a nonsteroidal anti-inflammatory drug (NSAID) used for treating pain, fever, and inflammation[1] fever (including postimmunisation fever), mild-to-moderate pain (including pain relief after surgery), painful menstruation, osteoarthritis, dental pain, headaches and pain from kidney stones. It is used for inflammatory diseases such as juvenile idiopathic arthritis and rheumatoid arthritis.[2][3] .It is also used for pericarditis and patent ductus arteriosus[4].

The pKa values of ligands and stability constants of the complexes with some hydroxamic acids a comparative study of three different potentiometric methods was reported by Senthilnithy [5]. The data obtained by pH-metric method were analyzed by three standard methods namely, Bjerrum's method, Irving and Rossotti method, and Sarkar and Kruck method. Acetohydroxamic acid, $CH_3CONHOH$, forms highly stable complexes with vanadium (V) and vanadium (IV) in 1: 1, 1: 2 and 1 : 3 mole ratios [6]. The stability of these complexes can be determined in terms of thermodynamic parameters; ΔG , ΔH and ΔS . The preliminary data, obtained through pH titration at various temperatures, was processed. Another study on the stability constant of the transition metal complexes with some medicinally important compounds was reported by Chaudhari [7]. The formation of bioligand complexes of some medicinal drugs with Co (II), Ni (II) and Cu (II) ions were investigated. The formation const. of Ranitidine Hydrochloride and 6-methoxy naphthaldehyde has been carried out pH metrically in aq. soln. at 30°C, at 0.1 M fixed ionic strength [7]. The method of Calvin and Bjerrum [8,9] as adopted by Irving and Rossotti [10] has been employed to determine log K value

Material and methods:

The titration were carried out in a 100 ml Pyrex glass beaker kept in a water bath maintain at constant temperature. Chemicals and ligand used were of analytical grade. Ligand solutions were prepared in ethanol. Metal salt solutions were prepared by dissolving the corresponding metal salt in twice distilled deionized water and standardized by standard volumetric methods. The free hydrogen ion concentrations were measured with a combined glass electrode attached to digital pH meter model-361; the accuracy of pH meter was 0.01 at three temperatures (25 ± 0.1 , 30 ± 0.1 and $35 \pm 0.1^\circ \text{C}$) and at an ionic strength of 0.1 mol L^{-1} (KNO_3). The pH meter was calibrated with suitable buffers before use.

Result and discussion :

Calvin – Bjerrum Titration :

The experimental procedure involved pHmetric titration of,
 Free acid (0.01M) titration,
 Free acid (0.01 M) and ligand (0.05M) titration.
 Free acid (0.01M), ligand (0.05M) and metal ion (0.01M) against std. NaOH solution.

The ionic strength of all solutions were maintain constant (0.1M) by adding appropriate quantity of 1M KNO_3 solution.

The titration were carried out in a 100 ml Pyrex glass beaker kept in a water bath maintain at constant temperature ($25, 30$ and 35°C) nitrogen gas was purged for chemically inert atmosphere. The readings were recorded for each addition of 0.1ml. The graphs of volume of alkali added against pH were plotted.

The titration curve of the acid and the ligand deviates at about pH 3.0 and then increase up to pH 12.0. The deviation between acid curve from ligand curve for the systems showed the dissociation of H^+ from-COOH groups of the ligands.

Proton – ligand formation Number (n^{A}):

Proton – ligand formation number (n^{A}) were calculated by Irving and Rossotti expression.

$$n^{\text{A}} = \gamma - (E^0 + N) (V_2 + V_1) / (V^0 + V_1) T_L^0$$

Where, V^0 = Initial volume of solution (50 ml)

N = Normality of sodium hydroxide

T_L^0 = Concentration of ligand in 50 ml solution

E^0 = Initial concentration of free acid (HNO_3)

γ = Number of dissociable proton from ligand

V_1 and V_2 – Volume of alkali consumed by acid and ligand on same pH

Metal- ligand formation number (n):

The deviation of ($A + L + M$) curve from ($A + L$) started at about pH 3.5, It indicate the commencement of complex formation.

Metal – Ligand formation number (n) was calculated by following expression.

$$n = (E^0 + N) (V_3 - V_2) / (V^0 + V_2) (T_M^0) \times n^{\text{A}}$$

V^0 = Initial volume of solution (50 ml)

N = Normality of sodium hydroxide

T_M^0 = Concentration of the metal ions

n^{A} = Proton – ligand formation number

E^0 = Initial concentration of free acid (HNO_3)

Where, V_2 and V_3 – volume of NaOH consumed by ligand and metal ions at same pH.

Metal – Ligand Formation curves:

Formation Curves were plotted between n and p^H . The metal-ligand stability constants were determined by half integral method

Half Integral Method:-The metal- ligand stability constants (\log, k_1 values) are calculated from formation curves. The values of $n = 0.5$ which corresponds to value of Pk . The values of metal – ligand stability constants i.e. $\log k$ for all the systems were presented in tables.

Table: Stability constants and thermodynamic parameters of Zr^{4+} with Ibuprofen

System	Temperature ($^{\circ}C$)	pKa	logK	$-\Delta H$ ($KJmol^{-1}$) At $30^{\circ}C$	$-\Delta G$ ($KJmol^{-1}$)	$-\Delta S$ ($KJmol^{-1}deg^{-1}$) At $30^{\circ}C$
Zr(IV)Ibuprofen	25	5.7446	7.265	5.50	41.452	0.12
	30	5.7385	7.185		41.683	
	35	5.7215	7.170		42.280	

Conclusion:

The results obtained from the pH metric measurements, the values of pKa were found to decrease with increasing temperature. The values of the thermodynamic functions ΔG , ΔH and ΔS were calculated. The values of stability constants reveal that the stability constants decrease with increasing temperature, along with the pKa value.

Acknowledgment :

Author are thankful to Principal ShriR.R.Lahoti Science College, Morshi for providing necessary laboratory facility

References:

1. The American Society of Health-System Pharmacists. Retrieved Jan 2016.
2. Joint Formulary Committee (2013). British National Formulary (BNF) (65 ed.). London, UK: Pharmaceutical Press. pp. 665, 671. ISBN 978-0-85711-084-8.
3. Rossi, S, ed. (2013). Australian Medicines Handbook (2013 ed.). Adelaide: The Australian Medicines Handbook Unit Trust. ISBN 978-0-9805790-9-3.
4. "Ibuprofen". The American Society of Health-System Pharmacists. Retrieved 3 April 2011.
5. Senthilnithy R, de Costa M.D.P, and H.D. Gunawardhana. H.D.
6. J. Natural Sci. Foundation Sri Lanka 2008, 36: 1916. Ali K., N. Fatima N and Z.T. Maqsood Z.T, Scientia 2005 Iranica, 12, 311
7. Chaudhari V.T., Ubale M. Band. Farooqui M. Journal of Ultra Chemistry, 2009, 5: 219
8. Calvin M. and Wilson K. W J. of American Chemical Society., 2003(1945). 67
9. J. Bjerrum, On the tendency of metal ions toward complex formation. In 6th Meeting of Scandinavian chemists, Lund, Sweden (1947).
10. Irving H.M. and Rossotti H.S, J. of the Chemical Society, (1954) 2904.
